# Coenzym Q<sub>10</sub> HPLC Kit

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|                                                         |           |
|---------------------------------------------------------|-----------|
| <b>1. INTENDED USE</b>                                  | <b>14</b> |
| <b>2. SUMMARY AND EXPLANATION OF THE TEST</b>           | <b>14</b> |
| <b>3. PRINCIPLE OF THE TEST</b>                         | <b>16</b> |
| <b>4. MATERIAL SUPPLIED</b>                             | <b>17</b> |
| <b>5. MATERIAL REQUIRED BUT NOT SUPPLIED</b>            | <b>17</b> |
| <b>6. PREPARATION AND STORAGE OF REAGENTS</b>           | <b>18</b> |
| <b>7. PRECAUTIONS</b>                                   | <b>18</b> |
| <b>8. SPECIMEN COLLECTION AND PREPARATION</b>           | <b>18</b> |
| <b>9. ASSAY PROCEDURE</b>                               | <b>18</b> |
| PROCEDURAL NOTES                                        | 18        |
| SAMPLE AND STANDARD PREPARATION                         | 19        |
| CHROMATOGRAPHIC CONDITIONS                              | 20        |
| <b>10. TREATMENT OF THE COLUMN</b>                      | <b>20</b> |
| <b>11. RESULTS</b>                                      | <b>20</b> |
| CALCULATION                                             | 20        |
| TYPICAL CHROMATOGRAM                                    | 22        |
| <b>12. LIMITATIONS</b>                                  | <b>22</b> |
| <b>13. QUALITY CONTROL</b>                              | <b>22</b> |
| EXPECTED VALUES                                         | 22        |
| CONTROLS                                                | 23        |
| <b>14. PERFORMANCE CHARACTERISTICS</b>                  | <b>23</b> |
| PRECISION AND REPRODUCIBILITY                           | 23        |
| LINEARITY                                               | 23        |
| DETECTION LIMIT                                         | 23        |
| <b>15. DISPOSAL</b>                                     | <b>23</b> |
| <b>16. TROUBLESHOOTING</b>                              | <b>24</b> |
| <b>17. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b> | <b>25</b> |

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Manual


# Coenzyme Q<sub>10</sub> HPLC Kit

*For the determination of Coenzyme Q<sub>10</sub> (Ubiquinone) in  
EDTA-whole blood, serum and EDTA-plasma*

Valid from 19.02.2013


KC 1700

|         |
|---------|
| CAL     |
| INT STD |
| CTRL 1  |
| CTRL 2  |

Immundiagnostik AG

## 1. INTENDED USE

The *Immundiagnostik* Assay is intended for the quantitative determination of Coenzyme Q10 (Ubiquinone) in EDTA-whole blood, serum and EDTA-plasma. This Assay is designed for *in vitro* diagnostic use only.

## 2. SUMMARY AND EXPLANATION OF THE TEST

Ubiquinone was first isolated in the 50th from the group of Prof. Green (Wisconsin). The function was investigated by Prof. Mitchel, who received the Nobel price for his research on the oxydative phosphorylation pathway.

Ubiquinone is a coenzyme which is represented in every cell and the whole metabolism. It is built up by a chinonic ring and an isoprenic sidechain. In humans ubiquinone can be synthesized and taken up by nutrition.

Ubiquinone has two different pysiological functions.

- 1.) Component of the energy metabolism
- 2.) Radical scavanger

During the reduction of oxygen in the oxydative phosphorylation 3 Mol ATP are generated. In this reduction electrons are tranferred from NADPH to oxygen via 6 different redox systems. Ubiquinone is the less abundant redox system in the membrane of the mitochondria. Because of the low amount it is the speed controlling redox-system in the energy metabolism. Normally the amount of ubiquinone is sufficient, but with growing age and exposure to sunlight it is reduced to 50 %.

Ubiquinone has a high amount of carbon doublebonds and therefore a higher potential of reduction than vitamin C or vitamin E. Thereby it is the first line of defense against free radicals. Therefore ubiquinone is an optimal stabilizer of the ion channels of the membranes.

### Indications

- Determination of Ubiquinone status
- Cardiovascular disease
- Carcinogenesis

- Aging
- Burnout syndrome

### **3. PRINCIPLE OF THE TEST**

The first step in determining Coenzyme Q10 includes sample preparation. Therefore the internal standard is added to the calibrator, controls and samples. Higher molecular substances are removed by precipitation and centrifugation. The Coenzyme Q10 in the supernatant is then extracted by an organic solvent, which is evaporated afterwards. The residue is resuspended in ethanol p.A. 100 µl of the resuspended solution are injected into the HPLC system.

The separation via HPLC follows an isocratic method at 30 °C using a reversed phase column. One run lasts 15 minutes. The chromatograms are recorded by a UV detector. The quantification is performed with the delivered calibrator; the concentration is calculated via integration of the peak areas by the internal standard method.

The kit includes all reagents for preparation and separation of the samples, except the column.

#### **Summary**

Besides many other parameters the advantage of HPLC method lies in the simultaneous handling of many analytes in a single test. The HPLC system enables even laboratories without experience in high performance liquid chromatography to use this technique for clinical routine determination in a quick and precise manner. Unlike immuno assays with up to six calibrators per test, a one-point calibration is sufficient to calibrate the test system. It is possible to automate the sample application and calculation of the results so that even higher numbers of samples can be handled nearly without control.

#### 4. MATERIAL SUPPLIED

| Cat. No   | Content        | Kit Components                      | Quantity    |
|-----------|----------------|-------------------------------------|-------------|
| KC 1700LM | MOPHA          | Mobile phase                        | 1000 ml     |
| KC 1700KA | CAL            | Calibrator (lyophilized 0.5 ml)     | 5 vials     |
| KC 1700IS | INTSTD         | Internal standard                   | 110 ml      |
| KC 1700RE | RECSOL         | Reconstitution solution             | 5 ml        |
| KC 1700VL | DIL            | Dilution solution                   | 85 ml       |
| KC 1700EX | EXTSOL         | Extraction solution                 | 220 ml      |
| KC 1700ET | ETHANOL        | Ethanol p.A.                        | 20 ml       |
| KC 1700KO | CTRL1<br>CTRL2 | Control 1 and 2, 250 µl lyophilized | 2 x 3 vials |

HPLC column (KC 1700 RP) as well as individual components can be ordered separately from Immundiagnostik. Please ask for the price list of the individual components.

#### 5. MATERIAL REQUIRED BUT NOT SUPPLIED

- 10 ml screw capped glass vials (Pyrex)
- Centrifuge
- Mixer (Vortex)
- Various pipettes
- HPLC with UV-detector
- Reversed phase C<sub>18</sub>-column
- Evaporation unit

## 6. PREPARATION AND STORAGE OF REAGENTS

- Reconstitute the **calibrator** (CAL) in **0.5 ml** reconstitution solution (RECSOL). Take aliquots and store at -20°C for at least 14 days. Avoid repeated freeze-thaw cycles. The coenzyme Q10 concentration might have minor variation from lot to lot.
- Reconstitute the **controls** (CTRL1, CTRL2) in **250 µl** reconstitution solution (RECSOL).
- Calibrator (CAL), internal standard (INTSTD) and controls (CTRL1, CTRL2) are stable at -20 °C up to the expiry date stated on the label.
- All other test reagents are stable at 2-8 °C up to the expiry date stated on the label.

## 7. PRECAUTIONS

- For *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Reagents should not be used beyond the expiration date shown on kit label.

## 8. SPECIMEN COLLECTION AND PREPARATION

EDTA-whole blood, serum or EDTA-plasma is suited for this test system.

The samples are stable for 1 day at 2-8°C in the dark. For longer storage samples, should be frozen at -20°C.

## 9. ASSAY PROCEDURE

### Procedural notes

- Quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may

influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.

- The assay should always be performed according the enclosed manual.
- Bring dilution solution (DIL) to room temperature before use to dissolve precipitations.

## Sample and standard preparation

Pipet into 10 ml screw capped glass vials

**200 µl** sample, calibrator (CAL) or control (CTRL1, CTRL2)

+

**800 µl** dilution solution (DIL)

Mix well (10 s vortex mixer)

Add **1 ml** Internal standard (INTSTD)

mix well (10 s vortex mixer)

Add **2 ml** extraction solution (EXTSOL)

Mix for **2 min** on a vortex mixer and centrifuge for **10 min** at 3.000 g.  
The upper layer is removed and evaporated.

Resuspend the residue in **150 µl** ethanol p.A. (ETHANOL)

Inject **100 µl** of the resuspended Solution into the HPLC system.

## Chromatographic conditions

|                          |                            |
|--------------------------|----------------------------|
| <b>Column material:</b>  | Bischof Prontosil AQ; 5 µm |
| <b>Column dimension:</b> | 125 mm x 4 mm              |
| <b>Flow rate:</b>        | 0.8 - 1.2 ml/min           |
| <b>UV-Detection :</b>    | 275 nm                     |
| <b>Temperature:</b>      | 30 °C                      |
| <b>Injection volume:</b> | 100 µl                     |
| <b>Running time:</b>     | 15 min                     |

It is recommended that a guard column is used to extend column life.

Use mobile phase (MOPHA) for autosampler and injection valve wash.

## 10. TREATMENT OF THE COLUMN

The column can be stored in the mobile phase (MOPHA) after the analysis.

Before use, the system should be equilibrated with 50 ml mobile phase (MOPHA): Run first 20 ml without column, and then the remaining 30 ml with integrated column.

## 11. RESULTS


### Calculation

$$\frac{\text{Peak height sample} \times \text{Concentration of the calibrator}}{\text{Peak height internal standard in the sample}} \times F = \text{Concentration sample}$$

$$F = \frac{\text{Peak height internal standard in the calibrator}}{\text{Peak height calibrator}}$$

## Typical chromatogram

## 12. LIMITATIONS

## 13. QUALITY CONTROL

### Expected values

**EDTA-whole blood: 0.67 – 0.99 µg/ml (mean ± 2 SD)**

**EDTA-plasma: 0.83 – 1.43 µg/ml (mean ± 2 SD)**

It is recommended that each laboratory should establish its own normal range. Above mentioned values are only for orientation and may vary from other published data.

## Controls

Control samples or serum pools should be analyzed with each run of calibrators and patient samples. Results generated from the analysis of the control samples should be evaluated for acceptability using appropriate statistical methods. In assays in which one or more of the quality control sample values lie outside the acceptable limits, the results for the patient sample may not be valid.

## 14. PERFORMANCE CHARACTERISTICS

### Precision and reproducibility

|                        |                    |          |
|------------------------|--------------------|----------|
| <b>Intra-Assay CV:</b> | 4.4 % (0.66 µg/ml) | [n = 12] |
|------------------------|--------------------|----------|

|                        |                    |          |
|------------------------|--------------------|----------|
| <b>Inter-Assay VK:</b> | 6,6 % (0,31 µg/ml) | [n = 10] |
|------------------------|--------------------|----------|

|  |                    |          |
|--|--------------------|----------|
|  | 4,5 % (0,89 µg/ml) | [n = 10] |
|--|--------------------|----------|

### Linearity

up to 10 µg/ml

### Detection limit

0.02 µg/ml

## 15. DISPOSAL

The mobile phase (MOPHA), extraction solution (EXTSOL), internal standard (INTSTD), and ethanol p.A. (ETHANOL) must be disposed as non-halogenated solvent.

Please refer to the appropriate national guidelines.

## 16. TROUBLESHOOTING

| Problem                  | Possible reason                                     | Solution                                                                                            |
|--------------------------|-----------------------------------------------------|-----------------------------------------------------------------------------------------------------|
| No signal                | No or defect connection to evaluation system        | Check signal cord and connection                                                                    |
|                          | Detector lamp is altered                            | Change lamp                                                                                         |
| No peaks                 | Injector is congested                               | Check Injector                                                                                      |
| Double peaks             | Dead volume in fittings and / or column             | Renew fittings and / or column                                                                      |
| Contaminating peaks      | Injector dirty                                      | Clean injector                                                                                      |
|                          | Contamination at the head of the column             | Change direction of the column and rinse for 30 min at low flow rate (0.2 ml/min) with mobile phase |
|                          | Air in the system                                   | Degas pump                                                                                          |
|                          | Autosampler vials contaminated                      | Use new vials or clean them with methanol                                                           |
| Broad peaks, tailing     | Precolumn / column exhausted                        | Use new precolumn / column                                                                          |
| Variable retention times | Drift in temperature                                | Use a column oven                                                                                   |
|                          | Pump delivers imprecise                             | Check pump, degas the system                                                                        |
|                          | System is not in steady state yet                   | Rinse system mobile phase for 15 min                                                                |
| Baseline is drifting     | Detector lamp did not reach working temperature yet | Wait                                                                                                |
|                          | Detector lamp is too old                            | Renew lamp                                                                                          |

|                        |                                   |                                      |
|------------------------|-----------------------------------|--------------------------------------|
|                        | System is not in steady state yet | Rinse system mobile phase for 15 min |
|                        | Pump delivers imprecise           | Check pump, degas the system         |
| Baseline is not smooth | Pump delivers imprecise           | Check pump, degas the system         |
|                        | Detector flow cell is dirty       | Clean flow cell                      |

## 17. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- The test components contain organic solvents. Contact with skin or mucous membranes has to be avoided.
- All reagents in the test package are for research use only.
- The reagents should not be used after the date of expiry stated on the label.
- Do not interchange different lot numbers of any kit component within the same assay.
- Quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.