# Homocystein HPLC Kit

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|                                                         |           |
|---------------------------------------------------------|-----------|
| <b>1. INTENDED USE</b>                                  | <b>14</b> |
| <b>2. SUMMARY AND EXPLANATION OF THE TEST</b>           | <b>14</b> |
| <b>3. PRINCIPLE OF THE TEST</b>                         | <b>14</b> |
| <b>4. MATERIAL SUPPLIED</b>                             | <b>15</b> |
| <b>5. MATERIAL REQUIRED BUT NOT SUPPLIED</b>            | <b>16</b> |
| <b>6. PREPARATION AND STORAGE OF REAGENTS</b>           | <b>16</b> |
| <b>7. PRECAUTIONS</b>                                   | <b>17</b> |
| <b>8. SPECIMEN COLLECTION AND PREPARATION</b>           | <b>17</b> |
| <b>9. ASSAY PROCEDURE</b>                               | <b>17</b> |
| PROCEDURAL NOTES                                        | 17        |
| SAMPLE AND STANDARD PREPARATION                         | 18        |
| CHROMATOGRAPHIC CONDITIONS                              | 18        |
| <b>10. TREATMENT OF THE COLUMN</b>                      | <b>19</b> |
| <b>11. RESULTS</b>                                      | <b>19</b> |
| CALCULATION                                             | 19        |
| TYPICAL CHROMATOGRAM                                    | 19        |
| <b>12. LIMITATIONS</b>                                  | <b>20</b> |
| <b>13. QUALITY CONTROL</b>                              | <b>20</b> |
| EXPECTED VALUES                                         | 20        |
| CONTROLS                                                | 20        |
| <b>14. PERFORMANCE CHARACTERISTICS</b>                  | <b>21</b> |
| PRECISION AND REPRODUCIBILITY                           | 21        |
| LINEARITY                                               | 21        |
| DETECTION LIMIT                                         | 21        |
| RECOVERY                                                | 21        |
| <b>15. DISPOSAL</b>                                     | <b>21</b> |
| <b>16. TROUBLESHOOTING</b>                              | <b>22</b> |
| <b>17. REFERENCES</b>                                   | <b>23</b> |
| <b>18. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b> | <b>23</b> |

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
# Homocysteine HPLC Kit

For the determination of Homocysteine in plasma and serum

Valid from 18.08.2010

KC 2801

|         |
|---------|
| CAL     |
| INT STD |
| CTRL 1  |
| CTRL 2  |

## **1. INTENDED USE**

The Immundiagnostik Assay is intended for the quantitative determination of homocysteine in plasma and serum. This Assay is designed for in vitro diagnostic use only.

## **2. SUMMARY AND EXPLANATION OF THE TEST**

Homocysteine exists in plasma as protein-bound and, so-called, free forms. About 70 percent of the total plasma homocysteine is bound to albumin. The free forms include homocystine (homocysteine disulphide) and mixed disulphide. Just minute amounts of reduced homocysteine are found in plasma. Total homocysteine in Plasma, including protein-bound and free forms, has been referred as homocyst(e)ine [H(e)]. Studies for peripheral vascular, cerebrovascular and coronary artery disease (CAD), evaluating H(e) as a cardiovascular risk factor, have been consistently demonstrated that there are more persons with high H(e) among patients than among healthy persons. Increased concentrations of free forms of homocysteine have also generally been reported among patients with occlusive artery disease, with one notable exception of a negative finding in relation to CAD.

## **3. PRINCIPLE OF THE TEST**

The first step in determining homocysteine includes sample preparation with additional derivatisation. Therefore the internal standard is added to the calibrator, controls and samples. Reduction and derivatisation is carried out in one step (10 min at 60 °C). During the reduction, homocysteine is cleaved from albumin and homocystine is reduced into two homocysteines. The derivatisation reagent transforms Homocysteine into a fluorescent product. Higher molecular substances are removed by precipitation and centrifugation. 20 ml of the supernatant are injected into the HPLC system.

The separation via HPLC follows an isocratic method at 30 °C using a reversed phase column. One run lasts 5 minutes. The chromatograms are recorded by a fluorescence detector. The quantification is performed with the delivered plasma calibrator; the concentration is calculated via integration of the peak areas by the internal standard method.

The application of Homocysteine for HPLC makes a fast and precise determination of the amino acid possible in an easy way. The kit includes all reagents for preparation and separation of the samples, except the column.

## Summary

Besides many other parameters the advantage of HPLC method lies in the simultaneous handling of many analytes in a single test. The HPLC system enables even laboratories without experience in high performance liquid chromatography to use this technique for clinical routine determination in a quick and precise manner. Unlike immuno assays with up to six calibrators per test, a one-point calibration is sufficient to calibrate the test system. It is possible to automate the sample application and calculation of the results so that even higher numbers of samples can be handled nearly without control.

## 4. MATERIAL SUPPLIED

| Cat. No   | Content          | Kit Components                           | Quantity    |
|-----------|------------------|------------------------------------------|-------------|
| KC 2801LM | MOPHA            | Mobile phase                             | 1 x 1000 ml |
| KC 2801KA | CAL              | Calibrator (lyophilized, 1 ml)           | 1 vial      |
| KC 2801IS | INT STD          | Internal standard (lyophilized, 12 ml)   | 1 vial      |
| KC 2801RE | RECSOL           | Reconstitution solution                  | 1 x 25 ml   |
| KC 2801RL | REDSOL           | Reduction solution (lyophilized, 2.4 ml) | 1 vial      |
| KC 2801DL | DER              | Derivatisation solution                  | 1 x 12 ml   |
| KC 2801FR | PREC             | Precipitation reagent                    | 1 x 12 ml   |
| KC 2801KO | CTRL 1<br>CTRL 2 | Control 1 and 2; 250 µl lyophilized      | 2 x 3 vials |

HPLC column (KC 2801RP) as well as individual components can be ordered separately from Immundiagnostik. Please ask for the price list of the individual components.

## 5. MATERIAL REQUIRED BUT NOT SUPPLIED

- 1.5 ml reaction tubes (Eppendorf)
- Centrifuge
- Ice water bath
- Various pipettes
- HPLC with Fluorescence-detector
- Reversed phase C<sub>18</sub>-column
- Water bath /Thermo shaker
- Vortex Mixer

## 6. PREPARATION AND STORAGE OF REAGENTS

- Reconstitute the **calibrator** (CAL) in **1 ml** reconstitution solution (RECSOL). Take aliquots and store them at -20°C. Avoid repeated freeze-thaw cycles. The concentration of homocysteine might have minor changes from lot to lot.
- Reconstitute the **controls** (CTRL1,2) in **250 µl** reconstitution solution (RECSOL).
- Reconstitute the **internal standard** (INT STD) in **12 ml** reconstitution solution (RECSOL). **Take aliquots** and store them at **-20°C**. Avoid repeated freeze-thaw circles. The concentration of homocysteine might have minor changes from lot to lot.
- Reconstitute the **reduction solution** (REDSOL) in **2.4 ml** reconstitution solution (RECSOL) and store at **2-8 °C**. The reconstituted reduction solution is stable for 3 months.
- All other test reagents are stable at **2-8 °C**, up to the date of expiry stated on the label.

## 7. PRECAUTIONS

- For in vitro diagnostic use only.
- This product contains human source material which was tested and found to be non-reactive to HBsAg, anti-HIV-1/2, and anti-HCV. Since no method can offer complete assurance that hepatitis B virus, HIV-1/2, HVC or other infectious agents are absent, these reagents should be handled as if potentially infectious.
- The precipitating reagent contains acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped out immediately with copious quantities of water. Do not breathe vapor and avoid inhalation.
- Reagents should not be used beyond the expiration date shown on kit label.

## 8. SPECIMEN COLLECTION AND PREPARATION

Since the concentration in serum samples is depending on clotting time, EDTA-plasma should be preferred as sample. Please store the sample directly after blood collection at 2-8°C. After blood collection, the sample should be centrifuged within the next 30 minutes (conditions: 2000 g, 10 min., 2-8°C). Afterwards the sample should be stored at 2-8°C (shipment is possible for 24 h at room temperature).

For long time storage (**> 1 week**), **keep the samples at -20 to -80°C.**

## 9. ASSAY PROCEDURE

### Procedural notes

- The quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure, that are not coordinated with the producer, may influence the test results. Immundiagnostik can therefore not be held reliable for any damage resulting from this.
- The assay should always be performed due to the manual which is given in the kit.

## Sample and standard preparation

**50 µl sample, CAL** (calibrator) or **CTRL 1, 2** (control 1 and 2)

+

**50 µl INT STD** (internal standard)

+

**20 µl REDSOL** (reduction solution)

+

**100 µl DER** (derivatisation reagent)

Incubate for **10 min** at 60 °C and afterwards stop the reaction in an ice water bath.

Add **100 µl** precipitating solution

Mix well, incubate for **10 min** at 2-8 °C and centrifuge for **10 min** at 10.000 g. The prepared sample is now stable for min. 6 days at 2-8 °C.

Inject **20 µl** of the supernatant into the HPLC system.

## Chromatographic conditions

|                                                                                            |                                                                          |            |        |
|--------------------------------------------------------------------------------------------|--------------------------------------------------------------------------|------------|--------|
| <b>Column material:</b>                                                                    | MZ Inertsil ODS-2; 5 µm                                                  |            |        |
|                                                                                            | MZ PerfectBond ODS-2; 5 µm                                               |            |        |
|                                                                                            | Bischoff Prontosil Eurobond; 5 µm                                        |            |        |
| <b>Column dimension:</b>                                                                   | 125 mm x 4 mm                                                            |            |        |
| <b>Flow rate:</b>                                                                          | 0.7 - 1.0 ml/min (please refer to the quality certificate of the column) |            |        |
| <b>Temperature:</b>                                                                        | 30 °C                                                                    |            |        |
| <b>Detection :</b>                                                                         | Fluorescence:                                                            | Excitation | 385 nm |
|                                                                                            |                                                                          | Emission   | 515 nm |
| <b>Injection volume:</b>                                                                   | 20 µl                                                                    |            |        |
| <b>Running time:</b>                                                                       | 5 min                                                                    |            |        |
| <b>Immundiagnostik recommends to use a guard-column to enlarge lifetime of the column.</b> |                                                                          |            |        |

## 10. TREATMENT OF THE COLUMN

After analysis, the column should be flushed with 30 ml aqua bidest. (1 ml/min) and stored in 50% methanol in aqua bidest (approx. 30 ml, flow 0.6 ml/min). Before use, the system should be equilibrated with ca. 30 ml eluent.


## 11. RESULTS

### Calculation

$$\frac{\text{Peak height sample} \times \text{Concentration of the calibrator}}{\text{Peak height internal standard in the sample}} \times F = \text{Concentration sample}$$

$$F = \frac{\text{Peak height internal standard in the calibrator}}{\text{Peak height calibrator}}$$

### Typical chromatogram

## 12. LIMITATIONS

Since the concentration in serum samples is depending on clotting time, EDTA-plasma should be preferred as sample.

Don't use whole blood.

## 13. QUALITY CONTROL

### Expected values

Homocysteine values: (Resch "Homocystein" ISBN 3-920328-17-5)

|        |        |                                            |
|--------|--------|--------------------------------------------|
| < 15   | μmol/l | Normal value                               |
| 15-30  | μmol/l | Lack of vitamins                           |
| 31-100 | μmol/l | Indication of heterocygote homocysteinemia |
| > 100  | μmol/l | Indication of homocygote homocysteinemia   |

After a methionine loading test an increase in homocysteine concentration is normal. Concentrations more than 30 μmol/l are indicating a heterocygote homocysteinemia.

We recommend each laboratory to establish its own normal range. The values mentioned above are only for orientation and can deviate from other published data.

### Controls

Control samples should be analyzed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

## 14. PERFORMANCE CHARACTERISTICS

### Precision and reproducibility

|                        |                     |          |
|------------------------|---------------------|----------|
| <b>Intra-Assay CV:</b> | 2.5 % (8.9 µmol/l)  | [n = 12] |
|                        | 1.6 % (16.5 µmol/l) | [n = 12] |
| <b>Inter-Assay CV:</b> | 5.2 % (8.7 µmol/l)  | [n = 12] |
|                        | 2.9 % (17.4 µmol/l) | [n = 12] |

### Linearity

up to 500 µmol/l

### Detection limit

0.6 µmol/l

### Recovery

98 %

## 15. DISPOSAL

The MOPHA (mobile phase), REDSOL (reduction solution), INT STD (internal standard), and DER (derivatisation solution) must be disposed as non-halogenated solvents. The PREC (precipitation solution) can be neutralized with NaOH to pH 7.0 and disposed as salt solution.

**Important:** Reaction will produce heat, be careful!

Please refer to the appropriate national guidelines.

## 16. TROUBLESHOOTING

| Problem                  | Possible reason                                     | Solution                                                                                            |
|--------------------------|-----------------------------------------------------|-----------------------------------------------------------------------------------------------------|
| No signal                | No or defect connection to evaluation system        | Check signal cord and connection                                                                    |
|                          | Detector lamp is altered                            | Change lamp                                                                                         |
| No peaks                 | Injector is congested                               | Check Injector                                                                                      |
| Double peaks             | Dead volume in fittings and / or column             | Renew fittings and / or column                                                                      |
| Contaminating peaks      | Injector dirty                                      | Clean injector                                                                                      |
|                          | Contamination at the head of the column             | Change direction of the column and rinse for 30 min at low flow rate (0.2 ml/min) with mobile phase |
|                          | Air in the system                                   | Degas pump                                                                                          |
|                          | Auto sampler vials contaminated                     | Use new vials or clean them with methanol                                                           |
| Broad peaks, tailing     | Precolumn / column exhausted                        | Use new precolumn / column                                                                          |
| Variable retention times | Drift in temperature                                | Use a column oven                                                                                   |
|                          | Pump delivers imprecise                             | Check pump, degas the system                                                                        |
|                          | System is not in steady state yet                   | Rinse system mobile phase for 15 min                                                                |
| Baseline is drifting     | Detector lamp did not reach working temperature yet | Wait                                                                                                |
|                          | Detector lamp is too old                            | Renew lamp                                                                                          |
|                          | System is not in steady state yet                   | Rinse system mobile phase for 15 min                                                                |
|                          | Pump delivers imprecise                             | Check pump, degas the system                                                                        |
| Baseline is not smooth   | Pump delivers imprecise                             | Check pump, degas the system                                                                        |
|                          | Detector flow cell is dirty                         | Clean flow cell                                                                                     |

## 17. REFERENCES

1. Role of Plasma Homocyst(e)ine in Arterial Occlusive Diseases. (1994) Clinical Chemistry, Vol. 40 No. 6, pp. 857-858.
2. Association between Plasma homocysteine concentrations and extracranial carotid-artery stenosis. Selhub, J. et al. (1995). The New England Journal of Medicine, Feb. 95, Vol. 332 No. 5, pp. 286-291.
3. Plasma Homocyst(e)ine as a Risk Factor for Early Familial Coronary Artery disease. Wu, L. et al. (1994). Clinical Chemistry, Vol. 40 No. 4, pp. 552-561.
4. Homocyst(e)ine and arterial occlusive diseases. Malinow, M. (1994). Journal of Intern Medicine, Vol. 236, pp. 603-617.
5. Hormone replacement therapy may reduce high serum homocysteine in postmenopausal women. van der Mooren, M. (1994). European Journal of Clinical Investigation, Vol. 24, pp. 733-736.
6. Homocysteine metabolism in pregnancies complicated by neural-tube defects. Mills, J. et al. (1995). The Lancet, Jan. 1995, Vol. 345, pp. 149-151.

## 18. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- The test components contain organic solvents. Avoid contact with skin or mucous membranes.
- For in vitro diagnostic use only.
- Reagents should not be used beyond the expiration date shown on the kit label.
- Do not interchange different lot numbers of any kit component within the same assay.
- The guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.