Hepcidin 25 LC-MS/MS Tuning Kit

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Hersteller

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Hepcidin 25 LC-MS/MS Tuning Kit



For determination of the optimal LC-MS/MS settings

EU: IVD / CE

US: Research Use Only. Not for use in diagnostic procedures.

Valid from 13.05.2014

REF **KM4001**

IVD **CE**

Immundiagnostik AG, Stubenwald-Allee 8a, D-64625 Bensheim

Tel.: +49 6251 70190-0

Fax: + 49 6251 849430

e.mail: info@immundiagnostik.com

www.immundiagnostik.com

Table of Contents

1. INTENDED USE	11
2. MATERIAL SUPPLIED	11
3. MATERIAL REQUIRED BUT NOT SUPPLIED	11
4. PREPARATION AND STORAGE OF REAGENTS	11
5. ASSAY PROCEDURE	12
<i>Tuning</i>	12
<i>Sensitivity test</i>	13
6. TREATMENT OF THE COLUMN	14
7. RESULTS	14
8. DISPOSAL	15
9. PRECAUTIONS	15
10. TECHNICAL HINTS	15
11. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	15

1. INTENDED USE

In preparation for the Immundiagnostik Hepcidin 25 LC-MS/MS application, an instrument tuning is necessary. The tuning kit is intended for determination of the optimal LC-MS/MS settings as well as for testing the sensitivity of the instrument.

The MS/MS parameters recommended by Immundiagnostik AG are only for orientation and must be optimized for each instrument.

2. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
KM4000LA	MOPHA A	Mobile Phase A	1000 ml
KM4000LB	MOPHA B	Mobile Phase A	1000 ml
KM4000AC	ACTSOL	Activating reagent	2.5 ml
KM4001TU	TUSOL	Tuning solution	5 ml
KM4001TS	TUSIS	Tuning solution, internal Standard (lyophilized)	2 ml
KM4001SO	SOL A	Solution A	25 ml

The UPLC column with inline filter (KM4000RP) can be ordered separately from Immundiagnostik.

3. MATERIAL REQUIRED BUT NOT SUPPLIED

- Vortex mixer
- 2 ml glass vials, suitable for LC-MS/MS
- 500 ml measuring cylinder, suitable for LC-MS/MS
- Different pipettors
- LC-MS/MS instrument
- RP-C18 column, e.g. XSelect CSH C18 (2,1 x 50 mm), 2,5 µm

4. PREPARATION AND STORAGE OF REAGENTS

Mobile phases and solution A

Add 0.1 % activation reagent (ACTSOL) to the mobile phases (MOPHA A, MOPHA B) and solution A (SOL A) before use:

e.g. 500 ml MOPHA (MOPHA A or MOPHA B) + 500 µl ACTSOL

The obtained solutions are stable for 2 weeks.

Warning: The activation reagent (ACTSOL) must be added under the fume hood. All vials to be used must be absolutely clean and detergent-free as well as preferentially made of glass suitable for LC-MS/MS analysis.

Tuning solution internal standard

Dissolve tuning solution internal standard (TUSIS) in 2 ml of activated solution A. The reconstituted solution is stable for maximally 1 week.

Tuning solution

Tuning solution (TUSOL) can be stored at -20°C and used until the expiration date given on the label.

Other reagents

All other reagents are stable at 2–8°C up to the date of expiry stated on the label.

5. ASSAY PROCEDURE

Tuning

We recommend tuning with a T-fitting (UPLC/HPLC flow rate and a syringe pump). The tuning solutions (TUSOL, TUSIS) consist of a highly pure hepcidin 25 and an internal standard (human calcitonin gene related peptide), with a concentration of 1 µg/ml each.

HPLC/UPLC flow rate of activated MOPHA A: 0.2 ml/min

Flow rate of syringe pump with the tuning solutions (TUSOL, TUSIS): 10–20 µl/min

MS/MS method (e. g. for a Waters Quattro Premier XE Tandem mass spectrometer)

Mode:	MRM
Polarity:	ESI ⁺
Capillary (kV):	4
Cone (V):	var.
Extracor (V):	4
RF Lens (V):	0
Source temperature (°C):	130
Desolvation Temperatur (°C):	500
Cone Gas Flow (l/h):	300

Desolvation Gas Flow (l/h): 950
 Collision Gas Flow (ml/min): 0.25

MRM transitions (m/z)

Hepcidin-25 (Molecular weight: 2789.4 g/mol)

698.41 > 644.54 Cone voltage: 37 Collision energy: 23

698.41 > 354.09 Cone voltage: 37 Collision energy: 33

The detected precursor ion corresponds to the 4-fold charged mass.

Internal Standard (Calcitonin Gene Related Peptide human; Molecular weight: 3789,31 g/mol)

758.91 > 718.83 Cone voltage: 40 Collision energy: 25

758.91 > 689.5 Cone voltage: 40 Collision energy: 25

The detected precursor ion corresponds to the 5-fold charged mass.

Sensitivity test

After obtaining a stable precursor ion and the corresponding fragment ions, a dilution series of the tuning solution are prepared to test the linearity and sensitivity of the LC-MS/MS system.

Preparation of the dilution series

The tuning solution (TUSOL) is diluted with activated solution A (SOL A) as follows:

1.) Preparation of the highest concentration:

TUSOL 1:2 in SOL A

e.g. 500 µl TUSOL in 500 µl SOL A (500 ng/ml)

2.) Preparation of the dilution series:

			Concentration [ng/ml]
1000 µl	500 ng/ml		500
500 µl	500 ng/ml	500 µl SOL A	250
400 µl	250 ng/ml	600 µl SOL A	100
100 µl	100 ng/ml	900 µl SOL A	10
100 µl	10 ng/ml	900 µl SOL A	1
		1000 µl SOL A	Blank

50 µl of the samples can be injected directly in the LC-MS/MS system. A duplicate injection of the samples is recommended. For determination of the retention time, solely the highest concentration of 500 ng/ml can be injected in the first step.

Chromatographic conditions:

Column material: e.g. XSelect CSH C18 (2,1 x 50 mm), 2,5 µm

Column dimension: 2.1 x 50 mm

Flow: 0.4 ml/min

Column temperature: 35 °C

Inject volume: 50 µl

Run time: 6 min

Gradient:

Time	% A	% B
0 min	90	10
1,5 min	90	10
3,0 min	5	95
4,0 min	5	95
5,0 min	90	10
6,0 min	90	10

It is recommended that a guard column or inline filter is used to extend column life.

6. TREATMENT OF THE COLUMN

After the analysis, the column should be washed with ~20 ml of 50 % methanol. The column can be stored in 50 % methanol.

7. RESULTS

The obtained areas are plotted against the concentrations and a linear regression curve is generated. The lowest concentration of 1 ng/ml should be found to obtain the necessary sensitivity.

When the obtained results are satisfying, the hepcidin 25 LC-MS/MS kit of Immundiagnostik AG can be used for measurement of hepcidin 25 in serum.

8. DISPOSAL

The mobile phases (MOPHA A, MOPHA B), tuning solutions (TUSOL, TUSIS), activating reagent (ACTSOL) and solution A (SOL A) must be disposed as non-halogenated solvents.

Please refer to the appropriate national guidelines.

9. PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- The mobile phases (MOPHA A and MOPHA B) as well as tuning solutions (TUSOL, TUSIS) and solution (SOL A) contain organic solvents and must be handled with gloves, eye protection, and appropriate protective clothing in a hood.
- The ACTSOL (activation reagent) consists of a strong acid and must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.

10. TECHNICAL HINTS

- Do not mix different lot numbers of any kit component.
- Reagents should not be used beyond the expiration date shown on the kit label.
- The assay should always be performed according the enclosed manual.

11. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- All reagents in the kit package are for *in vitro* diagnostic use only.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.

- Warranty claims and complaints in respect of deficiencies must be lodged within 14 days after receipt of the product. The product shall be send to Im-mundiagnostik AG together with a written complaint.

Used symbols:

	Temperature limitation		Catalogue Number
	In Vitro Diagnostic Medical Device		Contains sufficient for <n> tests
	Manufacturer		Use by
	Lot number		